



MEDICUR® pro

Use for:

- pain relief

and depending on
the frequency:

- relaxing
- balancing
- activating

Our mobile pain therapist

The **MEDICUR®pro** offers **three different programs/frequencies** for **alleviation of pain**. Use it as required with its **relaxing, balancing or activating** function.

A randomized, placebo-controlled, double-blind study conducted with the **MEDICUR®** (the previous version of the **MEDICUR®pro**) at King's College Hospital in London could prove the **efficacy of 7.8 Hz for knee arthrosis** according to academic medicine criteria. The biometric evaluation demonstrated a **highly significant effect** and pain relief. If you are interested, you can download or request the study (www.magnetotherapy.de / info@ams-ag.de).

Besides the two different programs/frequencies of **3.0 and 20 Hz** you can select also **7.8 Hz**, the principle value of the **Schumann frequency** spectrum. The **MEDICUR®pro** always also generates the **Geomagnetic frequencies** and the **Solar frequencies** - regardless wich program you choose. You find more information about Schumann, Geomagnetic and Solar frequencies inside this brochure.



MEDICUR®pro...

pain soothing →

is an **easy-to-use magnetic field therapy device with three programs**. It has been used very successfully for many years to **relieve pain** and is approved throughout the EU (Medical Product Laws, Class II b).

supports the regulation →

activates the body's self-healing processes in a natural way and supports its self-regulation. As a possible result, you might need fewer drugs and thus can avoid their side effects.

handy →

is **handy** and hardly larger than a pack of paper tissues.

automatic shutdown →

is **simply placed onto or next to the body**, on the solar plexus or directly on the affected area. It **turns off automatically after 10 minutes**. If necessary, it can be switched on again immediately. Overdosing is not possible. It has a physiologically sphere of 80 cm / 31.5 inches.



MEDICUR® pro



*analgesic,
relaxing, balancing,
activating*

Generally, each of the preset programs provides a possible **pain relief**.

In addition, the different frequencies have proven their value for certain situations, e.g.:

Frequency	Examples of use
3,0 Hz	<p>relaxation</p> <ul style="list-style-type: none"> • benefits relaxation and sleep • supports in case of nervousness and sleeping troubles
7,8 Hz	<p>balancing / stabilizing</p> <ul style="list-style-type: none"> • gives power, supports calmness, balance and concentration • 7,8 Hz is the principal value of the Schumann frequencies and our most important „well-being frequency“ • can be used to energize any food containing water by putting the device next to it for around one minute
20 Hz	<p>energizing</p> <ul style="list-style-type: none"> • sets free power and energy • useful for exhaustion, tiredness, debility, lack of energy and more

Producer: Snowden Healthcare
The Green, Ruddington
Nottingham, NG11 6HH, UK

Technical data:

Magnetic field frequencies: 3.0 / 7.8 / 20 Hz

Physiological range: 40 cm / Ø 80 cm
(15.7 inches / Ø 31.5 inches)

Battery: 9 Volt block battery

Weight: 110 g (3.9 ounces)

Size: 11 x 7.6 x 2.5 cm
(4.3 x 3.0 x 1.0 inches)

Scope of delivery: Device with battery and one spare battery, Velcro® strap, manual

Warranty: 10 years

Order number: 104

To see the current USD price and place your order please check our website and online shop:
www.magnetotherapy.de



The devices that are developed by **Dr. W. Ludwig®** are identified with the **arrow symbol**. This symbol stands for generating an electromagnetic biofield as in undisturbed nature. The frequency spectrum of the 64 essential trace elements (Geomagnetic frequencies) is generated by a worldwide unique process.

Buy the MEDICUR® pro without any risk!



After purchasing you can return the device within **90 days with a money-back guarantee**. You will see this also on your invoice. If you return the device in perfect condition within the 90 days, we will give you a **full refund** (less shipping costs).

AMS GmbH

Advanced Medical Systems

Hauptstr. 26 • D-97990 Weikersheim

GERMANY

Tel. +49 (0)79 34 / 99 34 89-0

info@ams-ag.de • www.magnetotherapy.de

Please order...

by telephone: +49 (0) 7934 / 99 34 89 – 0
or in our online shop
www.magnetotherapy.de

Thank you